

*A complete, fully programmable blot processing station*



## The Freedom Rocker™ BlotBot® Automated Blot Processor

### KEY BENEFITS

#### ► IMPROVE CONSISTENCY

Microprocessor controls ensure precise timing and dispensing. Get the same great results every time you use it.

#### ► RESULTS IN HALF THE TIME

Blots are ready to review when you arrive at work. Start your next test before lunch. Shave months off your research project.

#### ► USE EXISTING PROTOCOLS

No need to change the way you process your blots. Programming is easy with a color graphic display.

#### ► CONSERVE ANTIBODY

Use only small volumes per blot. Recover antibody stock for reuse.

**“I’m one of the first users of the [original] Freedom Rocker at Amgen in Cambridge, MA. The machine gets LOTS of use from all of us here.”**

DR. BARBARA GRUBINSKA, AMGEN

**NEXT ADVANCE**  
*Innovative Lab Products for the Life Sciences*

Incubate and wash your blots unattended overnight or during a busy day with the BlotBot®, the next generation Freedom Rocker.

### Stop babysitting your Westerns!

This fully programmable blot processing station automates the tedious incubation and washing steps to process blots and gels. You can be confident that all steps will be performed with consistency, from one day to another and from person to person. All blots receive reagent or wash buffer at the same time in equal volumes.

Simply dispense individual reagents into each tray and recover them separately. The innovative curved incubation tray reduces the amount of reagent required. The Freedom Rocker BlotBot is perfect for conserving precious antibodies.


### User-friendly Interface

The control panel has a large graphic display and an intuitive layout, making it very easy to program your own protocols. No need to compromise on your optimized conditions, either. Programs can include different rocking speeds, number of washes, and times for each step.

### Proven Track Record

Freedom Rockers are built to last. In fact, many of the original Freedom Rockers from 2005 are still in use today.

## Processing a Western Blot



The Freedom Rocker frees up your time, automating dozens of liquid handling steps and ensuring consistency.

## Features


- Automated delivery and removal of buffers/reagents
- Dispense different antibodies into each compartment
- Specially designed trays require very little reagent per blot
- Recover valuable reagents separately
- Just plug into an electrical outlet and use on your lab bench
- No air connections or vacuum pumps required
- Drip proof quick-release connections
- Patented technology
- 2 year warranty
- 30-day money back guarantee

**“The Freedom Rocker is perfect for silver staining of gels for mass spectrometry. It has made an unbelievable impact on data generation for me in cancer protein screening.”**

DR. STEVE GRIFFITHS, ATLANTIC CANCER RESEARCH INSTITUTE

**“The blot washer you sent us is perfect! I love it!”**

RENEE CLARK, CANADIAN FOOD INSPECTION AGENCY, LETHBRIDGE AB



## For More Information, Contact

Next Advance, Inc.

Tel. 1.800.738.1681

info@nextadvance.com

www.nextadvance.com